

AXS API




AGENDA

1	Introduction	All
2	AXS API check in – feedback on the AXS API	All
3	Straumann AXS Ecosystem	Samuel
4	Suggested Vatech – Straumann AXS Integration	Mahesh
5	Next steps	All

AXS ECOSYSTEM

AXS IS ENTRY POINT TO ALL STRAUMANN DIGITAL ECOSYSTEM



Synchronisation of patient information and patient files



FEEDBACK ON AXS API

AXS Public Developer Portal Docs API's

APIs

- Authentication
- Orders
- Organization
- Patients**
- Patient Documents (4)
- Webhook System

Patients v1 OAS 3.0

The Straumann AXS® Patients enables you to create and manage the patient journey.

Straumann

Servers

<https://ppr.api.us.aks.straumann.com - Pre-production - US region> [Authorize](#)

POST /v1/patients Create a new Patient

This endpoint allows authorized users to create new patient records in AXS.

Parameters

Name	Description
Organization-ID * required	ID of the customer organization for which the data belongs to.

string(sauid)
(header)

Request body required [Try it out](#)

application/json

Example Value | Schema

```
{
  "firstName": "John",
  "middleName": "Alexander",
  "lastName": "Doe",
  "birthDate": "1990-08-15",
  "customPatientId": "P1744804721",
  "gender": "MALE",
  "dataConsent": true
}
```

Sez GET Get Pa POST search Klon PR POST cor POST Cre Cotreat Test Pro POST sul

HTTP AXS API V1 > Patients > Create Patient

POST [{{base_url}}/v1/patients](#)

Docs Params Authorization Headers 11 Body Scripts Settings

☐ none ☐ form-data ☐ x-www-form-urlencoded ☒ raw ☐ binary ☐ GraphQL JSON

```
1 {
2   "customPatientId": "uniqueidentiifer",
3   "firstName": "Sam",
4   "lastName": "Demo",
5   "middleName": "",
6   "birthDate": "",
7   "gender": "MALE",
8   "issuingAuthority": "API",
9   "dataConsent": true
10 }
```

Response History

COLLECTIONS

- AXS API V1
 - OAuth
 - Lab Order workflow
 - Patients
 - POST Create Patient**
 - GET Get Patient By Id
 - POST Search Patients
 - POST Create Document
 - POST FSS File Upload using Create Document Response
 - DEL delete Documents By Id
 - POST rename Document By Id
 - POST search Documents By Criteria
 - GET Download Patient File
 - webhook
 - Company Consent
 - POST consent request
 - GET consent unlink

SUGGESTED VATECH – STRAUMANN AXS INTEGRATION

CONNECTING ON PREMISE DEVICES TO CLOUD API



• How the Flow Works:

1. Vatech Desktop app/server connects only to Vatech backend (**Vatech Auth Server**)
 2. Vatech auth server securely stores:
 - client_id
 - client_secret
 3. Vatech auth server:
 - Fetches access tokens from AXS platform
 - Calls APIs on behalf of the desktop app.
- **Desktop app never sees or handles sensitive credentials**

Key Design Principles

- Desktop applications are public clients
- No secrets (client_id / client_secret) must be stored locally.
- All secure communication must go through a trusted backend.

Benefits of This Approach

- **Strong Security:** Secrets never exposed to end-user devices
- **Centralized Secret Rotation:** Rotate credentials without redeploying desktop apps
- **Reduced Risk Surface:** Prevents reverse engineering exploits